

When "My data only" is selected the system will search through the user's specific investigations and data.

(This allows the user to search their "my data" list. This will be useful if a scientist is involved in thousands of investigations.)

This allows the user to view what page they are currently on, and how many pages there are in total.

The box containing "1" will be an input field. This allows the user to jump to a page they want to. They can navigate between pages by using the 'Prev' and 'Next' text buttons.

The aim of the design is to make it clear for the users where they are and where they can go next.

Note: Pages will be ordered by date and displaying the most recent first. This instruction will change to notify the user that the output is ordered by date.

Informs the user of how many items are on the current page, including how many items are remaining.

Currently logged into ISIS. Facility login

Catalog search

Investigation name: Start date:

Instrument name: End date:

Run range: Keywords:

Advanced search My data only

Search results: n investigations found

Double click on an investigation to view related datafiles.

Investigation title	Instrument	Facility	Run range	Start date	End date
Super mad experiment	SANS	ISIS	14553-14556	11/01/1066 00:00	23/01/1066 00:00
Na8 310K 002_offset	SANS	ISIS	14558-14569	11/01/1066 00:00	23/01/1066 00:00
dhemb57.7a bob room T	SANS	ISIS	2040-2099	21/03/1066 00:00	31/01/1066 00:00
Na8 310K 002_offset	SANS	ISIS	2090-3000	21/03/1066 00:00	31/01/1066 00:00

< Prev | Page of N | Next > Displaying 500 of 15,642

Datafile information

Title: Experiment that proves something awesome. Instrument: SANS2D
Proposal: CAL_SANS2D_2012-09-29T13:40:42 Run range: 147,148,190

<input type="checkbox"/>	Datafile name	Filesize	Start date	End date
<input checked="" type="checkbox"/>	SXD01677.RAW	10MB	11/01/1999 00:11	21/01/1999 00:00
<input checked="" type="checkbox"/>	SXD01677.LOG	1MB	11/01/1999 10:24	21/01/1999 10:25
<input type="checkbox"/>	SXD01697.LOG	1MB	21/03/1996 00:00	26/03/1996 18:24
<input type="checkbox"/>	SXD01697.RAW	24MB	21/03/1996 00:00	26/03/1996 20:21
<input checked="" type="checkbox"/>	BOB02931.NXS	10GB	19/05/2010 19:20	28/09/2012 19:20